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Check the instruments for contamination and clean, if required.
Clean diamond instruments with the cleaning stone G9920.

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
Check the instruments for contamination and clean, if required.
Clean tungsten carbide instruments with the wire brush P6820.

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Check polishers for damage, concentricity and contamination.
Clean the polisher surface and restore concentricity using the dressing instrument P4060.

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GB User Information for rotating dental and surgical instruments

Carefully observing the safety recommendations ensures optimum working conditions and extends life of your dental and surgical instruments.

- Use adequate supply of water spray to avoid damages to natural teeth and fillings. This should be achieved by means of a turbine with three coolant spray nozzles which distribute the coolant over the entire length of the instrument head. A minimum coolant supply of 50 ml/min is necessary in the >1,500 min⁻¹ speed range in order to ensure reliable cooling of the material being worked on.
- Moistening of the surfaces to be prepared must not be impaired by an incorrect suction extraction system, as deflection of the coolant jet can also cause heat related damage.
- Additional cooling is required for FG instruments with a total length exceeding 19 mm and for instruments with an instrument head diameter over 2 mm, since the cooling spray fixtures are adapted to the standard values for total length and instrument head diameter.
- Immediately remove any instruments that are damaged (e.g. blank spots in diamond coating), bent or that no longer run concentrically.
- The recommended applied contact pressure range is 0.3 - 2 N (30 - 200 p). Excessive contact pressure must be avoided. Particularly when using FG instruments, dosed contact pressure (non-forced application) is of essential importance for the grinding procedure and its consequences.
- Insert the instruments carefully and without using force as far as they can go into the chuck of the turbine or right angle. Miniature instruments should be used in special mini turbines only; if you use them in normal turbine, insert them no further than the end of the cylindrical part.
- The turbine and right angle must be in flawless technical condition. Continuously check proper operating the nozzles and other outlets.
- Avoid canting or levering when grinding.
- Let the instruments reach operational speed before applying them to the surface of the material you are working on.
- Observe the recommended speeds indicated in the table page 11.11.

Disinfection, cleaning and sterilization of rotating dental and surgical instruments

- Rotating instruments have to be disinfected, cleaned, dried and sterilized before the first application on the patient and immediately after every use.
- Immediately place used instruments in a disinfecting and cleaning solution. Limitations and material compatibility must be observed.
- Please use only approved disinfectants suitable for use with rotary instruments. The disinfectant must comply with legal requirements of and be licensed in the land in question. Alcohol is not suitable for use and does not provide sufficient disinfection. Please comply with the manufacturer's regulations for procedures before and after disinfecting.
- If instruments are kept in the cleaning bath too long, their code color may come off. Strong chemical agents such as hydrochloric acid, hydrogen peroxide, etc. may attack the instruments and are not suited for cleaning.
- Then follow the instructions of the manufacturer of the disinfecting and cleaning solution. Check that there is no residue. If necessary, clean the instruments again. Treat steel drills with anti-corrosion spray.
- Prior to sterilization, place the instruments in stands that are suitable for sterilization.
- Sterilize the instruments.
- After sterilization keep the instruments protected against recontamination.
- **Indication** : Sterilisation for polishers at 134° for 5 min is recommended. You must not use alcohol!
- Brushes are single-use products .

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