NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
e-mail: nti@nti.de • http://www.nti.de

DE RPR-Wurzelstift-Aufbausystem

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GB RPR Root Post Restoration Systems

Recommendation for cleaning, disinfection and sterilization

Except for the plastic posts (distancing posts), all parts of the RPR Root Post System can be disinfected and sterilized. Root posts and instruments can be cleaned with current cleaning solutions and then sterilized. Before sterilization, clean posts and instruments under running water and dry with air pressure.

Sterilization temperature 134°C; exposure time 5 min.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 ← www.nti.de

Plastic posts must not be sterilized!

Plastic Material is not resistant to solvents and must not be disinfected with organic substances such as alcohol.

User information for RPR Ceramicor Posts (Casting instructions):

The melting range for Ceramicor cast on alloy is 1,400°C - 1,490°C. The RPR Ceramicor Posts have a diameter of 0.5 mm (yellow) to 1.1 mm (blue) on the spruing point. The casting temperature of the alloy, which will be casted on the Ceramicor Post must not exceed 1,400°C as to avoid the melting or deforming of the RPR Ceramicor post during casting.

RPR Root Canal Posts should be used after sterilization (description above) as follows:

- ① Endodontic treatment.
- ② Prepare a concave surface on the root face with a round diamond bur.
- ③ Select a post which ends approximately 2 mm before the root apex.
- ④ Use the correct size of root facer to drill 1 - 2 mm into the canal, do not lever, risk of breakage. Use a diamond bur (e.g. 835 008, -010) to cut a notch in the edge of the flat area. The lug helps prevent rotation.

- ⑤ Fit the correct depth gauge for the length of post to the root canal drill.
- ⑥ Enlarge the root canal by hand or with a handpiece.
Speed: 500 - 1,000 min⁻¹.
- ⑦ Roughen the root canal walls with the roughening instrument.
- ⑧ Clean and dry the root canal with paper points.
- ⑨ Coat the post with small particle, low viscosity cement - applying less to the tip - and press the post into place.
- ⑩ Allow to set and remove the excess cement.
- ⑪ Build up the core with composite or glass ionomer cement.
- ⑫ Take the impression and fabricate the restoration in the laboratory.

RPR Ceramicor and plastic posts for building up directly in the mouth:

When using cast on Ceramicor posts, adhere strictly to the casting instructions.

1. Refer to sections 1 - 8.
2. There are two methods for building up the core directly:

2.1 In the surgery (Application by Dentist):

Apply silicone or Vaseline to the cavity and insert the plastic or Ceramicor post. Build up the core in the mouth using a resin approved for oral use. Remove the post and core and have it cast in the laboratory.

2.2 In the laboratory:

Preparatory measures carried out by the dentist: Insert the plastic or Ceramicor post into the root canal and take the impression over it.

3. Once the post and core has been cast, cement it into place.

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