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**Reprocessing recommendation as per DIN EN ISO 17664
for rotary dental and surgical instruments**

General principles

All instruments must be cleaned, disinfected and sterilised before each use; this applies particularly to the first use after delivery because all instruments are supplied non-sterile (cleaning and disinfection after removal of the protective transport packaging; sterilisation after wrapping).

(Only process the ProxyStrip system after separating it into its constituent parts (handle, ProxyStrip)).

Efficient cleaning and disinfection is an essential prerequisite for effective sterilisation.

As part of your responsibility for the sterility of the instruments during use, as a matter of principle please ensure that only adequately validated methods specific for equipment and product are used for the cleaning/disinfection and sterilisation, that the equipment used (ultrasonic bath, steriliser) are regularly maintained and checked and that the validated parameters are complied with for every cycle.

Please follow as a matter of principle the legal regulations applicable in your country as well as the hygiene regulations in the medical practice or hospital. This applies particularly to the different specifications for effective prion inactivation.

Due to the intended purpose and the geometry/surface properties, we recommend classification as semi-critical B or critical B (depending on the type and location of the actual use).

Cleaning and disinfection

Principles

An automatic method (WD – washer-disinfector) should be used where possible for cleaning and disinfection. Because of the clearly lower effectiveness and reproducibility, a manual method – using an ultrasonic bath – should only be used if an automatic method is not available.

Pre-treatment

Immediately after use (within 2 h at most) coarse contamination must be removed from the devices.

Separate the ProxyStrip system into its constituent parts (handle, ProxyStrip) for this purpose.

To do so, use running water or a disinfectant solution; the disinfectant should be free of aldehyde (blood contamination may otherwise be fixed), have a tested effectiveness (e.g. DGHM/VAH or FDA/EPA approval or CE marking), be suitable for instrument disinfection and be compatible with the instruments (see 'Material resistance' section).

For manual removal of contamination only use a soft brush that you use solely for this purpose and never use metal brushes or steel wool.

For profiled instruments in particular, ensure that you have covered all surfaces and that there are no visible residues remaining and then rinse thoroughly (at least 1 min) under running water.

Instruments that have contamination adhering to their surfaces that is not soluble in water and cannot be manually removed must be sorted out at this point.

When cleaning internally cooled instruments, particular care must be taken. Thoroughly flush internally cooled instruments with fully deionised (DI) water. If the cooling holes are not open, replace the instruments.

Please note that the disinfectant used for the pre-treatment only serves for personal protection and does not replace the subsequent disinfection step to be carried out after completed cleaning.

Automatic cleaning/disinfection (WD)

When selecting the WD, ensure

- that the WD has as a matter of principle tested effectiveness (e.g. DGHM or FDA/EPA approval or CE marking as per DIN EN ISO 15883),
- that if possible a tested program for thermal disinfection (A_0 value ≥ 3000 or – for older devices – at least 5 min at 90°C) is used (for chemical disinfection, there is a risk of disinfectant residues on the instruments),
- that the program used is suitable for the instruments and includes sufficient rinsing cycles,
- that only sterile or low-microbe (max. 10 microbes/mL) and low-endotoxin (max. 0.25 endotoxin units/mL) water (e.g. purified water / highly purified water) is used for final rinsing,
- that the air used for drying is filtered (oil-free, low-microbe and low-particle) and
- that the WD is regularly maintained and checked.

When selecting the cleaning agent system to be used, ensure

- that it is as a matter of principle suitable for cleaning instruments made of metals and plastics,
- that – if thermal disinfection is not used – a suitable disinfectant with tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval or CE marking) is also used and that this is compatible with the cleaning agent used and
- that the chemicals used are compatible with the instruments (see 'Material resistance' section).

**Reprocessing recommendation as per DIN EN ISO 17664
 for rotary dental and surgical instruments**

The concentrations, temperatures and exposure times stated by the manufacturer of the cleaning agent and, if applicable, disinfectant as well as the guidelines for final rinsing must be complied with.

- Procedure:
1. Place the disassembled instruments into the WD using a small parts basket (if a basket of the right size is not available, a small parts basket does not necessarily have to be used for the handle of the ProxyStrip system). Ensure when doing so that the instruments do not make contact with one another.
 2. Start the program.
 3. Remove the instruments from the WD once the program has ended.
 4. Inspect and pack the instruments as soon as possible after they are removed from the disinfectant (see sections on 'Inspection', 'Maintenance' and 'Packing'), where necessary after further drying in a clean location.

Manual cleaning and disinfection

When selecting the cleaning agent and disinfectant to be used, ensure

- that they are in principle suitable for cleaning and disinfecting instruments,
- that the cleaning agent and disinfectant is suitable for ultrasound cleaning (no foam formation),
- that a cleaning agent and disinfectant with tested effectiveness (e.g., DGHM/VAH or FDA/EPA approval or CE marking) is used and
- that the chemicals used are compatible with the instruments (see 'Material resistance' section).

The prerequisite for using a combined cleaning agent/disinfectant is a very low initial load (no visible contamination) as the result of an effective preliminary cleaning of the instruments.

The concentrations, temperatures and exposure times stated by the manufacturer of the cleaning agent and disinfectant as well as the guidelines for final rinsing must be complied with. Only use freshly prepared solutions, only use sterile or low-microbe (max. 10 microbes/mL) and low-endotoxin (max. 0.25 endotoxin units/mL) water (e.g. purified water / highly purified water) and only use medical-grade compressed air for drying (oil-free, low-microbe and low-particle).

- Procedure:
- Cleaning
1. Disassemble the instruments as far as possible.
 2. Place the disassembled instruments into the cleaning bath (ultrasonic bath) so that the instruments are completely covered.
 Ensure when doing so that the instruments do not make contact with one another.
 3. Then activate the ultrasound for the specified exposure time (but not less than 5 min).
 4. Then remove the instruments from the cleaning bath and rinse them thoroughly at least three times (at least 1 min) with water.
 5. Dry instruments with medical-grade compressed air (with reprocessable compressed air gun, air pressure: max 10 bar, room temperature).
 6. Inspect the instruments (see sections on 'Inspection' and 'Maintenance').
 7. Pack the instruments as soon as possible after they are removed from the disinfectant (see section on 'Packing'), where necessary after further drying in a clean location.

Inspection

Inspect all instruments after cleaning or cleaning/disinfection (using a magnifying glass with 10× magnification) for corrosion, damaged surfaces, chips, deformations (e.g. bent and non-concentric instruments) as well as contamination and separate out any damaged instruments (for restrictions on the number of times an instrument can be reused, see 'Reusability' section). Instruments that are still contaminated must be cleaned and disinfected again.

Maintenance

Maintenance is not required.

Instrument oils must not be used.

**Reprocessing recommendation as per DIN EN ISO 17664
for rotary dental and surgical instruments**

Wrapping

(Only wrap the ProxyStrip system after separating it into its constituent parts (handle, ProxyStrip)).

Wrap the instruments in single-use sterilisation wrapping (single or double layer wrapping) that satisfies the following requirements:

- compliance with DIN EN 868-2ff/DIN EN ISO/ANSI AAMI ISO 11607
- suitability for steam sterilisation (temperature resistance to at least 138°C (280°F), sufficient steam permeability)
- adequate protection of the instruments or sterilisation wrapping against mechanical damage
- adequate size of the wrapping that ensures the seal is not placed under tension

Sterilisation

For the sterilisation only the following list of sterilisation methods may be used; other sterilisation methods are not permitted.

Steam sterilisation

- fractionated vacuum method at least **3-phase** (at least 4-phase for internally cooled instruments) with sufficient product drying. (The drying time depends primarily on aspects that are the sole responsibility of the user and must therefore be validated by the user.)
- steam steriliser that complies with DIN EN 13060 **type B** or DIN EN 285 or ANSI AAMI ST79
- validated in accordance with DIN EN ISO/ANSI AAMI ISO 17665 (valid commissioning and product-specific performance assessment)
- maximum sterilisation temperature 134°C (273°F; plus tolerance according to DIN EN ISO/ANSI AAMI ISO 17665)
- **sterilisation time (exposure time at the sterilisation temperature) at least 5 min at 132°C (270°F)/134°C**

| - Country | Fractionated vacuum method | Gravitation method |
|-----------------|---|---------------------------|
| Germany | at least 5 min at 134°C | not recommended/permitted |
| Other countries | at least 5 min at 132°C (270°F) / 134°C (273°F) | not recommended/permitted |

The flash sterilisation method is not permitted as a matter of principle.

Do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation, no chemiclaves and no plasma sterilisation.

Storage

After sterilisation the instruments must be stored in the sterilisation wrapping (must not be under tension) in clean, dry and dust-free conditions and protected against moisture and recontamination (room temperature approx. 20°C and 40–60% air humidity). The instruments must also be protected against chemicals, acids, sunlight, heat and extreme temperature fluctuations. Please note the shelf life that results from the validation of the sterilisation wrapping.

Material resistance

Ensure when selecting the cleaning agent and disinfectant that they do not contain the following ingredients:

- organic, mineral and oxidising acids
- strong alkalis (pH > 10.5 not permitted, only neutral or weakly alkaline cleaning agents are recommended)
- for polishing and ceramic instruments do not use alkaline cleaning agents or alcohol-based cleaning agents and disinfectants
- alcohols, ethers and ketones, benzine
- oxidants

General instructions:

- FlexiPol polishing discs and Profilax ECO polishers are single-use products and may only be manually disinfected before the first use. Reprocessing in a WD is not possible. FlexiPol polishing discs and Profilax ECO polishers cannot be sterilised.
- Reprocess polishers and ceramic abrasives in the instrument bath because the vibrations in the ultrasonic bath can be absorbed by the material.
- Steel burs made of tool steel are not suitable for either automatic reprocessing or steam sterilisation and can only be manually disinfected. The user should switch to appropriate hard metal instruments.
- Colour-anodised aluminium parts lose their colour when standard cleaning methods or WDs are used. During the reprocessing, specially coordinated cleaning agents and disinfectants must be used.
- Never clean any instruments and sterilisation trays with metal brushes or steel wool. Instruments and sterilisation trays must not be exposed to temperatures greater than 138°C (280°F)!
- Root posts made of glass fibre are single-use products. They can be cleaned with alcohol. They cannot be sterilised or disinfected in a WD. The use of autoclaves is not recommended. After cleaning do not touch with bare hands.

**Reprocessing recommendation as per DIN EN ISO 17664
for rotary dental and surgical instruments**

Reusability

The following values are empirical values for the reusability (product service life) of the product groups listed below:

| | | | |
|---------------------------------|-------|-------------------------------|-----------|
| Diamond instruments: | - 10× | Polishers: | - 10× |
| Instruments of stainless steel: | - 10× | Hard metal instruments: | - 15× |
| | | Endo instruments: wide canals | - max. 6× |
| | | moderate canals | - max. 3× |
| | | narrow canals only use | 1× |

The end of the product service life (time at which the reprocessing can no longer be considered safe) is defined by faulty instruments e.g. with missing diamond grit, blunt/broken blades, fractured working parts, corroded surfaces, bent instruments, etc. The end of the product service life is determined by the operator who then sorts out the faulty instruments.

Single-use items (labelled with ) are not permitted to be reused. Safe use cannot be guaranteed when reusing these products because there is a risk of infection and/or the safety of the products can no longer be assured.

Using damaged or contaminated instruments is the responsibility of the user.
Any liability is excluded with failure to follow these instructions.

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