NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 nti@nti.de • www.nti.de

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GB Instructions for use: Dental diamond grinding instruments

Content:

1. Scope of application
2. Intended use
3. Indication
4. Contraindication
5. Intended patient group
6. Intended users
7. User environment
8. Cleaning, disinfection, sterilisation
9. Reusability
10. Storage
11. Correct use
12. Warnings
13. Disposal
14. General instructions
15. Adverse effects
16. Conduct in case of serious incidents
17. Colour coding
18. Description of the symbols used
19. Warranty
20. Returns

Please read the instructions for use carefully before using the instruments and follow them at all times. NTI-Kahla GmbH shall not be liable for any damages due to non-compliance with these instructions or as a result of use that is incompatible with the system.

1. Scope of application:

These instructions for use apply to Dental diamond grinding instruments for the preparation of crown stumps, the preparation of cavities, work on fillings, root planing. These are Class IIa medical devices. They are marked with .

2. Intended use:

All Dental diamond grinding instruments covered by these instructions for use are supplied non-sterile, are intended for multiple use, require reprocessing for this purpose and are designed for use with corresponding dental motors (active medical devices). The Dental diamond grinding instruments are divided into the following four areas of application.

2.1. Diamond grinding instruments for crown stump preparation

The instruments are intended for intraoral work and/or the removal of tooth structure and/or dental materials within the context of the preparation techniques for crown and bridge restoration and are marked as follows:

2.2. Diamond grinding instruments for cavity preparation

The instruments are intended for intraoral work and/or the removal of tooth structure and/or dental materials within the context of the preparation techniques for a cavity and are marked as follows:

2.3. Diamond grinding instruments for work on fillings

The instruments are intended for intraoral work with dental technical and/or dental restorative materials (final fit) and are marked as follows:

2.4. Diamond grinding instruments for root planing

The instruments are intended for removing subgingival calculus and are marked as follows:

3. Indications:

3.1. Indications for diamond grinding instruments for crown stump preparation

Indications include dental caries, disorders of tooth development, other diseases of dental hard structures, faulty occlusion, processing dental materials, as well as other dental diseases.

3.2. Indications for diamond grinding instruments for cavity preparation

Indications include dental caries, other diseases of dental hard structures, processing dental materials, as well as other dental diseases.

3.3. Indications for diamond grinding instruments for work on fillings

Indications include faulty occlusion, processing dental materials, as well as other diseases of teeth and the periodontium.

3.4. Indications for diamond grinding instruments for root planing

Indications include diseases of the periapical tissue, gingivitis as well as other diseases of the periodontium.

4. Contraindications:

4.1. Contraindications for the diamond grinding instruments for crown stump preparation as well as for cavity preparation

Contraindications include cleaning of teeth within the context of prophylaxis, (maxillofacial) surgery, diseases of the periapical tissue, gingivitis and diseases of the periodontium.

4.2. Contraindications of the diamond grinding instruments for work on fillings

Contraindications include cleaning of teeth within the context of prophylaxis, (maxillofacial) surgery, diseases of the periapical tissue, gingivitis and diseases of the periodontium as well as diseases of the dental hard structures.

4.3. Contraindications for diamond grinding instruments for root planing

Contraindications include cleaning of teeth within the context of prophylaxis, (maxillofacial) surgery, dental caries, disorders of tooth development, other diseases of dental hard structures, faulty occlusion as well as processing of dental materials.

5. Intended patient groups:

The intended patient groups include children with permanent dentition, pregnant and lactating women, and adults of all ages.

6. Intended users:

Dental diamond grinding instruments are intended for use in dentistry and may only be used by persons who are familiar with the use of these instruments on the basis of their dental training and experience.

7. User environment:

Dental diamond grinding instruments are used at dental treatment stations employing appropriate dental motors.

8. Cleaning, disinfection, sterilisation:

Dental diamond grinding instruments are not supplied sterile. Prior to first use and immediately after use, the Dental diamond grinding instruments must be reprocessed according to the reprocessing recommendations according to EN ISO 17664. This can be accessed at www.nti.de. Upon request, the recommendations can be provided in a printed format within seven days.

9. Reusability:

Reusability depends on wear and possible damage incurred during use. The instruments must be inspected as part of reprocessing. Sort out worn and damaged instruments and do not reuse them.

10. Storage:

The Dental diamond grinding instruments must be stored in a dry place, protected from direct sunlight, at room temperature (approx. 20 °C/68 °F), in their original packaging or protected from (re-)contamination. Observe possible restrictions of the sterile barrier systems used.

11. Correct use:

To remove the instruments from the packaging, peel off the blister paper. Forcing the instruments through the blister paper damages the instruments.

Prior to initial use and immediately after use, the instruments must be cleaned, disinfected, dried, checked and sterilised.

After cleaning or cleaning/disinfection, the instruments must be checked for corrosion, surface damage, chipping, deformation as well as contamination using a 10x magnifying glass. Contaminated instruments must be cleaned and disinfected again. Damaged instruments are to be discarded and must not be re-used.

The use of personal protective equipment is necessary for all applications.

The Dental diamond grinding instruments are intended exclusively for use in motors conforming to EN ISO 14457. They must be in perfect technical condition. Clamp the instrument as deeply as possible into the drive system without using force and check that it is firmly seated. Secure the instrument against dropping, swallowing and aspiration during intraoral use.

Never exceed the maximum speed specific for the instrument. Start the instrument before coming into contact with the material to be machined.

Ensure sufficient cooling of at least 50 ml/min.

Wetting/cooling must not be impaired.

For instruments with a total length of >19 mm and/or a working part diameter >1.8 mm, additional external cooling is required.

A maximum contact pressure of 2 N must not be exceeded during application.

12. Warnings:

The Dental diamond grinding instruments are very small parts with pointed geometries, sharp edges and cutting blades. Swallowing or aspirating instruments or parts of instruments can cause injuries to the respiratory or the digestive tracts. Secure the instruments against dropping, swallowing or aspiration when using them intraorally. If an instrument or part of an instrument are swallowed or aspirated - for example due to breakage of an instrument - provide first aid, seek medical attention immediately and follow Section „16. Conduct in case of serious incidents“.

Needle-prick injuries can occur when inserting and removing instruments into or from the motor's collet.

During use, the instruments must not be jammed due to excessive contact force, tilted or be used as a lever. This results in damage to the dental hard structure, to dental materials and the instruments, as well as build-up of heat, overheating, chipping of the working parts, damage to the pulp as well as an increased risk of breakage of the instruments. Breakage of instruments can lead to swallowing or aspiration.

Insufficient cooling as well as too high speeds can lead to irreversible damage of the pulp or the surrounding tissue.

Exceeding the instrument-specific maximum speed can cause the instruments to break. Breakage of instruments can lead to swallowing or aspiration.

To remove the Dental diamond grinding instruments from the packaging, the blister paper must be peeled off. Forcing the instruments out may result in damage to the instruments and make them unusable.

Non-compliance with the instructions for correct use (see Section 11) as well as the warnings (see Section 12) can lead to premature wear, destruction of the instruments or the motor, damage to tissues and hazards to the person using the instrument, the patient as well as third parties.

13. Disposal:

The Dental diamond grinding instruments fall under the waste legislation category within the EU of "pointed and sharp medical instruments", the so-called "sharps". They may possibly be bacterially contaminated and must be handled with particular care to prevent needle-stick injuries of the users and third parties. The instruments must be disposed of in waste containers that cannot be penetrated and are unbreakable and impervious to mois-

ture and which do not release their contents under pressure or impact or when dropped. The waste containers must be clearly identifiable by their colour and labelling  (biological risks). The waste containers must be disposed of in accordance with the laws of the respective country by specialist companies. The European Waste Code is 180101.

14. General instructions:

Further instructions for use as well as the recommendations for reprocessing can be found at www.nti.de.

Upon request, the recommendations can be provided in a printed format within seven days.

15. Adverse effects:

Due to the pointed geometries, sharp edges and cutting edges, swallowing or aspirating instruments or parts of instruments can cause injuries to the respiratory or the digestive tracts.

When using the Dental diamond grinding instruments, unintentional damage or injury to tissue structures such as teeth, mucous membranes and bones as well as necrotisation of same can occur.

Incorrect reprocessing of the Dental diamond grinding instruments can lead to inflammatory reactions and infections in patients.

Furthermore, handling of the instruments can cause injuries to the user due to the pointed geometries and sharp edges.

Therefore, please observe these instructions for use at all times - in particular Section „12. Warnings - as well as the reprocessing recommendation“

16. Conduct in case of serious incidents:



In case of serious incidents, please proceed according to the legal provisions and report them to the competent authority as well as to NTI-Kahla GmbH.

17. Colour coding:

The grit sizes of the Dental diamond grinding instruments are indicated by the coloured ring on the shank:

supercoarse grit – **SC** – black ring
 coarse grit – **C** – green ring
 medium grit – **M** – blue ring or no marking
 fine grit – **F** – red ring
 superfine grit – **SF** – yellow ring
 ultrafine grit – **UF** – white ring

18. Description of the symbols used:

 Angle	 Catalogue number
 Spray cooling / use with water spray	 ISO number
 Ultrasonic bath	 Lot number / enables tracking of the corresponding production batch
 Washer-disinfector for thermal disinfection	 Manufacturer
 Sterilisable in steam steriliser (autoclave) at the specified temperature	 Date of manufacture
 Medical device	 Biological risks
 Attention	 Keep away from sunlight
 Observe electronic instructions for use (eIFU)!	 Keep dry
 Maximum permissible speed	 Non-sterile
 Package unit	 Open blister

19. Warranty:

NTI-Kahla GmbH excludes all warranty claims and accepts no liability whatsoever for direct or consequential damage resulting from misuse, improper use, improper preparation as well as non-compliance with the instructions for use.

20. Returns:

Any products for return may only be sent back to NTI-Kahla GmbH after disinfection/sterilisation has been performed and made clearly visible as such.

NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 nti@nti.de • www.nti.de

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