(DE) (GB) (FR)  
 (IT) (ES) (RU)

NTI-Kahla GmbH • Rotary Dental Instruments  
 Im Camisch 3 • D-07768 Kahla/Germany  
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29  
 e-mail: nti@nti.de • http://www.nti.de

## **(DE) Produktinformation KR Tisch**

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## **(GB) KR Table Product Information**

Dear customer,

Please read the product information carefully before using the KR Table. The guarantee does not cover damage caused by incorrect or inappropriate use.

### **Contents:**

- 1. Applications of the system**
- 2. Recommended use of the KR Table**
- 3. Recommended use with polishers**
  - Dressing polishers
  - CeraGlaze porcelain polisher
- 4. Technique for checking rotary instruments**
- 5. Technique for locks**

#### **Check that the delivery is complete and includes the correct number of items:**

1 KR Table

1 Dressing diamond, medium grit (no colour ring)

1 Dressing diamond, coarse grit (green ring)

1 Allen key for securing the dressing diamond

Transport and storage box (wooden)

1 Product information leaflet

#### **1. Applications of the system:**

- Dressing polishers to a precise angle of 0°, 1°, 2°, 4° or 6°.
- Fabricating locks.

#### **2. Recommended use of the KR Table:**

Regular care and maintenance increases the durability and service life of the table. Clean the table with compressed air or a steam cleaner, particularly in the holes for the dressing diamonds. Remove ingrained dirt with a small brush and washing-up liquid or in a boiling-out unit.

Check the relevant channel in the table and the shank for dirt before inserting the dressing diamond. Clean the channels with a bottle brush. Finger tighten the locking screw with the Allen key supplied (Order no. KR-1003).

Check the base of the KR Table for dirt and fix it in the required position with the magnet of the milling unit.

The dressing diamonds are subject to a certain degree of wear and tear, which varies depending on the polisher used. The round shoulder is particularly affected. Check the dressing diamond for damage before each use. The shank of the dressing diamond has grooves, which allow it to be used in four positions. This makes it possible to use any side of the dressing diamond.

Dressing diamonds can be reordered using the following order numbers:

Dressing diamond, medium, Order no. KR-2001M

Dressing diamond, coarse, Order no. KR-2002C

The dressing diamonds are precision instruments. To maintain accuracy and avoid damage, the instruments should only be used in the KR Table. A diamond dressing stone (Order no. P4060) is available for dressing and cleaning polishers manually.

#### **3. Recommended use with polishers:**

The polishers should be checked before each milling procedure:

- Check whether the required angle is still accurate.
- The surface of the polisher should be smooth without grooves, as these could detrimentally affect the accuracy. We recommend dressing the polisher again beforehand with the medium dressing diamond.
- The surface of the polisher must be clean. If there is any dirt, clean the surface with the medium dressing diamond.

#### **Dressing all types of polishers:**

Insert the polisher in the milling unit. Set the motor speed to 5000 – 8000 rpm. Move along the coarse dressing stone horizontally for dressing to an angle of 0°

applying light pressure.

Vertical dressing (from top to bottom) is used to attain the correct angle of 1°, 2°, 4° or 6°. The angle is produced with the coarse (green ring) dressing diamond. Once the polisher has been dressed to the correct angle, the rounded end is produced by dressing horizontally.

Then smooth the surface of the polisher (except the green CeraGlaze, see below) with the medium dressing diamond. The polisher is now ready to use.

#### **PLEASE NOTE:**

If excessive pressure is applied when dressing the polisher, it is possible that the tip will deflect and not be at the correct angle. This can have a detrimental effect on results, particularly with softer polishers.

#### **CeraGlaze porcelain polisher from NTI-Kahla GmbH:**

The grit size of the green CeraGlaze porcelain polisher is coarser than the diamond coating of the medium dressing diamond. The green CeraGlaze should not therefore be used on the medium diamond dressing stone.

The CeraGlaze polisher has a very stable bond. The high stability of the bond facilitates milling and extends the service life. The hardness produces greater frictional heat. A cosmetic sponge, water spray or a syringe should also be used for cooling when polishing porcelain primary units. Do not focus polishing on one point but keep the polisher moving continuously.

#### **4. Technique for checking rotary instruments**

The KR Table can be used for checking the concentricity and angle of various instruments. Insert the dressing diamond with the medium grit size in the required position in the KR Table and lock it in place. Switch off the drive motor of the spindle and move the instrument manually along the dressing diamond. Rotate the spindle with one finger to check the instrument.

#### **5. Technique for locks**

There are five axes with diameters of 1.0 mm, 1.2 mm, 1.3 mm, 1.4 mm and 1.5 mm on the KR Table. These diameters correspond to the most common axis diameters of customized swinglocks.

Place the lock on the relevant axis. Hold the lock in position with one finger so that the milling surface is above the indented area in the KR Table. Mill the edge using a parallel rotary cutter to a precise angle of 90°.

The individual axes are used in the holes in the indentations.

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