


Frästechnik/Milling Technique / Technique du fraisage 2,35 mm/3 mm

Wachsfräser konisch rund
Wax spiral cutter, conical round
Fraise spirale à cire conique à bout arrondi

Laboratório • Laboratorio • Laboratório

HF356WS (200)



Größe/Size/Taille \varnothing 1/10 mm	023	031	040
Handstück kurz • Handpiece short • pièce à main courte L mm	13,0	13,0	13,0
Winkel/Angle/Angle	2°	4°	6°

Shank ISO REF 1

HP \varnothing 2,35 mm 500 103 200 364 ... HF356WS-...	023	031	040
HP \varnothing 3,00 mm 500 123 200 364 ... HF356WS-123-...	023		

Konusfräser Kreuzverzahnung grob, rund
Cone cutter coarse cross cut, round end
Fraise conique denture croisée grosse, bout arrondi

HF356XFR (200)



Größe/Size/Taille \varnothing 1/10 mm	023	031
Handstück kurz • Handpiece short • pièce à main courte L mm	13,0	13,0
Winkel/Angle/Angle	2°	4°

Shank ISO REF 1

HP \varnothing 2,35 mm 500 103 200 220... HF356XFR-...	023	031
HP \varnothing 3,00 mm 500 123 200 220... HF356XFR-123-...	023	

Konusfräser Einfachverzahnung rund
Cone cutter plain cut, round end
Fraise conique denture simple, bout arrondi

HF356SR (200)




Größe/Size/Taille \varnothing 1/10 mm	023
Handstück kurz • Handpiece short • pièce à main courte L mm	13,0
Winkel/Angle/Angle	2°

Shank ISO REF 1

HP \varnothing 2,35 mm 500 103 200 135... HF356SR-...	023
HP \varnothing 3,00 mm 500 123 200 135... HF356SR-123-...	023

Konusfräser Kreuzverzahnung rund
Cone cutter cross cut, round end
Fraise conique denture croisée, bout arrondi

HF356XR (200)



Größe/Size/Taille \varnothing 1/10 mm	023	031	040
Handstück kurz • Handpiece short • pièce à main courte L mm	13,0	13,0	13,0
Winkel/Angle/Angle	2°	4°	6°

Shank ISO REF 1

HP \varnothing 2,35 mm 500 103 200 190... HF356XR-...	023	031	040
--	------------	------------	------------

Konusfräser Einfachverzahnung rund mit Fasenschliff
Cone cutter plain cut, round end with chamfer ground section
Fraise conique denture simple, bout arrondi avec chanfrein spécial

HF356SRF (200)



Größe/Size/Taille \varnothing 1/10 mm	023
Handstück kurz • Handpiece short • pièce à main courte L mm	13,0
Winkel/Angle/Angle	2°

Shank ISO REF 1

HP \varnothing 2,35 mm 500 103 200 103... HF356SRF-...	023
HP \varnothing 3,00 mm 500 123 200 103... HF356SRF-123-...	023