

## Frästechnik/Milling Technique / Technique du fraisage 2,35 mm/3 mm

**Parallel Wachsspiralfräser rund**  
**Parallel spiral wax cutter, round end**  
**Fraise spirale à cire parallèle à bout arrondi**

**HF364WS (137)**

Größe/Size/Taille ø 1/10 mm	010	015	023
Handstück kurz • Handpiece short • pièce à main courte L mm	8,0	10,0	15,0

Shank

HP ø 2,35 mm 500 103 137 364... <b>HF364WS-...</b>	<b>010</b>	<b>015</b>	<b>023</b>
HP ø 3,00 mm 500 123 137 364... <b>HF364WS-123-...</b>	<b>010</b>		

**Parallelfräser Kreuzverzahnung grob, rund**  
**Parallel cutter coarse cross cut, round end**  
**Fraise parallèle denture croisée grosse, bout arrondi**

**HF364XFR (137)**

Größe/Size/Taille ø 1/10 mm	010	015	023
Handstück kurz • Handpiece short • pièce à main courte L mm	8,0	10,0	15,0

Shank

HP ø 2,35 mm 500 103 137 220 ... <b>HF364XFR-...</b>	<b>010</b>	<b>015</b>	<b>023</b>
HP ø 3,00 mm 500 123 137 220 ... <b>HF364XFR-123-...</b>	<b>010</b>	<b>015</b>	<b>023</b>

**Parallelfräser Kreuzverzahnung rund**  
**Parallel cutter cross cut, round end**  
**Fraise parallèle denture croisée, bout arrondi**

**HF364XR (137)**

Größe/Size/Taille ø 1/10 mm	010	015	023
Handstück kurz • Handpiece short • pièce à main courte L mm	8,0	10,0	15,0

Shank

HP ø 2,35 mm 500 103 137 190 ... <b>HF364XR-...</b>	<b>010</b>	<b>015</b>	<b>023</b>
HP ø 3,00 mm 500 123 137 190 ... <b>HF364XR-123-...</b>	<b>010</b>	<b>015</b>	<b>023</b>

**Parallelfräser Einfachverzahnung rund, linksdrall**  
**Parallel cutter plain cut, round end, left twist**  
**Fraise parallèle denture simple, bout arrondi, à gauche**

**HF364RL (137)**

Größe/Size/Taille ø 1/10 mm	010	015	023
Handstück kurz • Handpiece short • pièce à main courte L mm	8,0	10,0	15,0

Shank

HP ø 2,35 mm 500 103 137 135 ... <b>HF364RL-...</b>	<b>010</b>	<b>015</b>	<b>023</b>
HP ø 3,00 mm 500 123 137 135 ... <b>HF364RL-123-...</b>	<b>010</b>	<b>015</b>	<b>023</b>


**Parallelfräser Einfachverzahnung rund, Linksdrall mit Fasenschliff**  
**Parallel cutter plain cut, round end, left twist, flat with chamfer ground section**  
**Fraise parallèle denture simple, bout arrondi, à gauche, avec chanfrein spécial**

**HF364RLF (137)**

Größe/Size/Taille ø 1/10 mm	010	015	023
Handstück kurz • Handpiece short • pièce à main courte L mm	8,0	10,0	15,0

Shank

HP ø 2,35 mm 500 103 137 103 ... <b>HF364RLF-...</b>	<b>010</b>	<b>015</b>	<b>023</b>
HP ø 3,00 mm 500 123 137 103 ... <b>HF364RLF-123-...</b>	<b>010</b>	<b>015</b>	<b>023</b>



HF364RLF-015