



Hartmetallfräser • Tungsten Carbide Cutters • Fraises en carbure de tungstène

CE Kreuzverzahnung standard • Standard Cross Cut • Denture croisée standard

Labor • Laboratory • Laboratoire

Größe/Size/Taille ø 1/10 mm	023	023	023	060	040	060	023	014	
L mm	14,0	16,0	16,0	12,5	7,8	12,2	16,0	3,5	
Shank	ISO	1							
HP	500 104 ...	201 190 ...	137 190 ...	292 190 ...	113 190 ...	263 190 ...	263 190 ...	116 190 ...	274 190 ...
REF HF ...		261CE-023	293CE-023	295CE-023	296CE-060	351CE-040	351CE-060	364CE-023	390CE-014



CN Einfachverzahnung standard • Standard Plain Cut • Denture simple standard





HF078CN-060

HF071CN-027

HF257CNR-060

Größe/Size/Taille ø 1/10 mm	027	031	040	050	060	060	031	040	045	
L mm	2,4	2,8	3,4	4,3	12,2	12,2	11,5	14,1	12,9	
Shank	ISO	1								
HP	500 104 ...	001 175 ...	001 175 ...	001 175 ...	001 175 ...	238 175 ...	257 175 ...	199 175 ...	201 175 ...	200 175 ...
REF HF ...		071CN-027	071CN-031	071CN-040	071CN-050	L077CN-060	078CN-060	079CN-031	079CN-040	079CN-045

Größe/Size/Taille ø 1/10 mm	060	009	010	040	060	060	060	060	
L mm	15,0	4,2	4,2	11,5	13,9	14,2	13,9	12,2	
Shank	ISO	5		1					
HP	500 104 ...	201 175 ...	196 175 ...	196 175 ...	274 175 ...	274 175 ...	275 175 ...	194 175 ...	263 175 ...
REF HF ...		079CN-060	138CN-009	138CN-010	251CN-040	251CN-060	L251CN-060	257CNR-060	351CN-060