


Hartmetallfräser • Tungsten Carbide Cutters • Fraises en carbure de tungstène

GE Kreuzverzahnung grob • Bulk Cross Cut • Denture croisée grosse

Labor • Laboratory • Laboratoire




Größe/Size/Taille ø 1/10 mm	045	060	060
L mm	14,1	13,9	12,2
Shank			
HP 500 104 ...	201 220 ...	274 220 ...	263 220 ...
REF HF ...	079GE-045	251GE-060	351GE-060



HF251 GE-060

QSCN Einfachverzahnung mit Querhieb • Coarse Straight Blade Cross Cut • Denture simple avec taille transversale




Größe/Size/Taille ø 1/10 mm	040	045	023	060	060	070
L mm	9,0	14,1	14,0	14,2	12,2	14,2
Shank						
HP 500 104 ...	237 176 ...	201 176 ...	201 176 ...	275 176 ...	263 176 ...	263 176 ...
REF HF ...	077QSCN-040	079QSCN-045	261QSCN-023	L251QSCN-060	351QSCN-060	351QSCN-070




HFL251 QSCN-060

DF Diamantschliff fein, angeraute Oberfläche
Diamond Fine Cut For Slightly Rough Surface
Denture diamantée fine à surface rugueuse




Größe/Size/Taille ø 1/10 mm	045	023
L mm	14,1	8,0
Shank		
HP 500 104 ...	201 141 ...	198 141 ...
REF HF ...	079DF-045	138DF-023

D Diamantschliff Prismenoberfläche
Diamond Cut With Prism Surface
Denture diamantée à surface prismatique



Größe/Size/Taille ø 1/10 mm	045
L mm	14,1
Shank	
HP 500 104 ...	201 191 ...
REF HF ...	079D-045

QF Spiralverzahnung fein • Spiral Fine Cut • Denture hélicoïdale fine



Größe/Size/Taille ø 1/10 mm	023	040	023	023	023	023	040	060	023	023
L mm	5,3	13,7	9,0	9,5	10,5	11,5	11,5	13,9	10,5	14,0
Shank										
HP 500 104 ...	237 134 ...	201 134 ...	141 134 ...	289 134 ...	199 134 ...	274 134 ...	274 134 ...	274 134 ...	185 134 ...	201 134 ...
REF HF ...	L077QF-023	079QF-040	129QF-023	139QF-023	199QF-023	251QF-040	251QF-060	251QF-060	257QF-023	261QF-023

Bitte Drehzahlbereiche beachten (Drehzahltable - Seite 11.15) / Please observe Recommended speeds (speed chart - page 11.15) / Respecter les vitesses de rotation recommandées (voir à la page 11.15)