




Hartmetallfräser • Tungsten Carbide Cutters • Fraises en carbure de tungstène

GE Kreuzverzahnung grob • Bulk Cross Cut • Denture croisée grosse

Labor • Laboratory • Laboratoire






Größe/Size/Taille ø 1/10 mm	045	060	060
L mm	14,1	13,9	12,2
Shank  			
HP 500 104 ...	201 220 ...	274 220 ...	263 220 ...
REF HF ...	079GE-045	251GE-060	351GE-060



HF251 GE-060

QSCN Einfachverzahnung mit Querhieb • Coarse Straight Blade Cross Cut • Denture simple avec taille transversale






Größe/Size/Taille ø 1/10 mm	040	045	023	060	060	070
L mm	9,0	14,1	14,0	14,2	12,2	14,2
Shank  						
HP 500 104 ...	237 176 ...	201 176 ...	201 176 ...	275 176 ...	263 176 ...	263 176 ...
REF HF ...	077QSCN-040	079QSCN-045	261QSCN-023	L251QSCN-060	351QSCN-060	351QSCN-070




HFL251 QSCN-060



DF Diamantschliff fein, angeraute Oberfläche
Diamond Fine Cut For Slightly Rough Surface
Denture diamantée fine à surface rugueuse




Größe/Size/Taille ø 1/10 mm	045	023
L mm	14,1	8,0
Shank  		
HP 500 104 ...	201 141 ...	198 141 ...
REF HF ...	079DF-045	138DF-023



D Diamantschliff Prismenoberfläche
Diamond Cut With Prism Surface
Denture diamantée à surface prismatique



Größe/Size/Taille ø 1/10 mm	045
L mm	14,1
Shank  	
HP 500 104 ...	201 191 ...
REF HF ...	079D-045

QF Spiralverzahnung fein • Spiral Fine Cut • Denture hélicoïdale fine








Größe/Size/Taille ø 1/10 mm	023	040	023	023	023	023	040	060	023	023
L mm	5,3	13,7	9,0	9,5	10,5	11,5	13,9	10,5	14,0	
Shank  										
HP 500 104 ...	237 134 ...	201 134 ...	141 134 ...	289 134 ...	199 134 ...	274 134 ...	274 134 ...	185 134 ...	201 134 ...	
REF HF ...	L077QF-023	079QF-040	129QF-023	139QF-023	199QF-023	251QF-040	251QF-060	257QF-023	261QF-023	

Bitte Drehzahlbereiche beachten (Drehzahltable - Seite 11.15) / Please observe Recommended speeds (speed chart - page 11.15) / Respecter les vitesses de rotation recommandées (voir à la page 11.15)

Hartmetallfräser • Tungsten Carbide Cutters • Fraises en carbure de tungstène



QF Spiralverzahnung fein • Spiral Fine Cut • Denture hélicoïdale fine








Größe/Size/Taille ø 1/10 mm	008	009	010
L mm	4,2	4,2	4,2
Shank  			
HP 500 104 ...	196 134 ...	196 134 ...	196 134 ...
 HF ...	138QF-008	138QF-009	138QF-010

Zum Feinausarbeiten von Fissuren
Optimal Cosmetic Contouring and Fine Works on Fissures
 Denture hélicoïdale fine pour la finition de sillons

QSF Spiralverzahnung superfein • Spiral Superfine Cut • Denture hélicoïdale superfine









Größe/Size/Taille ø 1/10 mm	040	023	008	016	023	023	060
L mm	13,7	8,0	4,2	8,0	8,0	8,0	14,2
Shank  							
HP 500 104 ...	201 137 ...	141 137 ...	196 137 ...	198 137 ...	198 137 ...	289 137 ...	275 137 ...
 HF ...	079QSF-040	129QSF-023	138QSF-008	138QSF-016	138QSF-023	139QSF-023	L251QSF-060





HF129QSF-023




GTi Titanverzahnung grob • GTi Special cutter for titanium • GTi Fraise spéciale pour le titane

Größe/Size/Taille ø 1/10 mm	040	023	016	023	023
L mm	13,7	8,0	8,0	8,0	8,0
Shank  					
HP 500 104 ...	201 194 ...	141 194 ...	198 194 ...	198 194 ...	289 194 ...
 HF ...	079GTi-040	129GTi-023	138GTi-016	138GTi-023	139GTi-023

FLE / LE Linkshänderfräser Spezialfräser, linksschneidend • L-Cut Special Cutters • Fraise spéciale pour gauchers

Größe/Size/Taille ø 1/10 mm	045	023	045	060
L mm	12,9	16,0	12,9	13,9
Shank  				
HP 500 104 ...	200 140 ...	292 140 ...	200 240 ...	274 240 ...
 HF ...	079FLE-045	295FLE-023	079LE-045	251LE-060

LE = KREUZVERZÄHNUNG GROB
 LE = CROSS CUT COARSE
 LE = DENTURE CROISÉE GROSSE

FLE = FEIN • FINE • FINE