


## Endodontie • Endodontics • Endodontie



### 178 ...

**Wurzelfüller ohne Sicherheitsspirale**  
**Root fillers without safety spiral**  
**Bourre-pâte sans spirale de sécurité**



Größe / Size / taille ø 1/100 mm

**025** **030** **035** **040**

Shank	ISO	REF				
RA L=21	340 204 671 458 ...	<b>178.21.204. ...</b>	4	<b>025</b>	<b>030</b>	<b>035</b> <b>040</b>
(1 piece each)	<b>Assortment</b>	<b>178.21.204.SO</b>	1	025	030	035 040
RA L=25	340 204 673 458 ...	<b>178.25.204. ...</b>	4	<b>025</b>	<b>030</b>	<b>035</b> <b>040</b>
(1 piece each)	<b>Assortment</b>	<b>178.25.204.SO</b>	1	025	030	035 040

⌚ opt. 800 rpm

Ohne Sicherheitsspirale. Winkelstück mit grüner Farbcodierung.  
 Without safety spiral. RA handpiece with green colour coding.  
 Sans spirale de sécurité. Contre - angle avec code couleur vert.

Endo Stopper  
 Endo Stops  
 Endo Stops





**REF ES 250**  250



### 178S ...

**Wurzelfüller mit Sicherheitsspirale**  
**Root fillers with safety spiral**  
**Bourre-pâte avec spirale de sécurité**



Größe / Size / taille ø 1/100 mm

**025** **030** **035** **040**

Shank	ISO	REF				
RA L=21	340 204 672 458 ...	<b>178S.21.204. ...</b>	4	<b>025</b>	<b>030</b>	<b>035</b> <b>040</b>
(1 piece each)	<b>Assortment</b>	<b>178S.21.204.SO</b>	1	025	030	035 040
RA L=25	340 204 673 458 ...	<b>178S.25.204. ...</b>	4	<b>025</b>	<b>030</b>	<b>035</b> <b>040</b>
(1 piece each)	<b>Assortment</b>	<b>178S.25.204.SO</b>	1	025	030	035 040

⌚ opt. 800 rpm

Winkelstück mit grüner Farbcodierung.  
 RA handpiece with green colour coding.  
 Contre-angle avec couleur vert.

