


Finierer • Finishing Carbide Burs • Instruments à finir

Praxis • Dental Surgery • Cabinet dentaire

H390 (274)

 Granate
 Grenade
 Grenade



Größe/Size/taile ø 1/10 mm	014	016	018
Länge-Length in mm - longueur	3,5	3,5	3,7

 Shank   

 FG  500 314 274 072... **H390 - ... -FG** **014(12)** **016(12)** **018(12)**

⚙️ opt. 160.000 rpm

H41 (001)

 Rund
 Round
 Rond


Größe/Size/taile ø 1/10 mm	014	016	018	023
Länge-Length in mm - longueur	1,2	1,4	1,6	2,06

 Shank   


US. No. 7004 7006 7008

 FG  500 314 001 071... **H41 - ... -FG** **014(12)** **016(12)** **018(18)** **023(20)**

 RA  500 204 001 071... **H41 - ... -RA** **018(18)** **023(20)**

⚙️ opt. FG 160.000 rpm. RA 30.000 rpm

H46 (254)

 Flamme
 Flame
 Flamme


Größe/Size/taile ø 1/10 mm	012	014	018	023
Länge-Length in mm - longueur	3,5	3,5	3,8	4,6

 Shank   

US. No. 7102 7104 7106 7108

 FG  500 314 254 072... **H46 - ... -FG** **012(12)** **014(12)** **018(12)**


 RA  500 204 254 072... **H46 - ... -RA** **018(12)** **023(12)**


⚙️ opt. FG 160.000 rpm, RA 30.000 rpm

H48L (249)

 Flamme lang
 Long Flame
 Flamme longue


Größe/Size/taile ø 1/10 mm	010	012
Länge-Length in mm - longueur	8,0	8,0

 Shank   

 FG  500 314 249 072... **H48L - ... -FG** **010(12)** **012(12)**


 FG  500 314 249 041... **H48LF - ... -FG** **012(20)**

 FG  500 314 249 031... **H48LUF - ... -FG** **012(30)**

⚙️ opt. 160.000 rpm

Keramik Finierer für die Laborturbine Ceramic Finishing Instruments for Laboratory Turbine Fraises à finir pour la céramique pour la turbine de laboratoire

H8503 (467)

 Keramik Finierer/ 3 - Kant
 Ceramic Finishing Bur/ 3-sided
 Fraise à finir/ 3 arêtes


Größe/Size/taile ø 1/10 mm	014
Länge-Length in mm - longueur	2,5


 Shank   

Arbeitsteil • Working part • Partie travaillante 9°

 FG 500 314 467 211... **H8503 - ... -FG** **014(3)**

⚙️ opt. 160.000 rpm

H8504 (467)

 ... / 4 - Kant
 ... / 4 - sided
 ... / 4 - arêtes


Größe/Size/taile ø 1/10 mm	012
Länge-Length in mm - longueur	2,5


 Shank   

Arbeitsteil • Working part • Partie travaillante 10°


 FG 500 314 467 212... **H8504 - ... -FG** **012(4)**

⚙️ opt. 160.000 rpm

H8506 (467)

 ... / 6 - Kant
 ... / 6 - sided
 ... / 6 - arêtes





Größe/Size/taile ø 1/10 mm	010
Länge-Length in mm - longueur	2,5

 Shank   

Arbeitsteil • Working part • Partie travaillante 12°

 FG 500 314 467 212... **H8506 - ... -FG** **010(6)**

 HP 500 104 467 212... **H8506 - ... -HP** **010(6)**
 ⚙️ opt. 160.000 rpm, HP 25.000 rpm

	...	Roter Ring/Red ring/Bague rouge	Standard	ISO No. 071 - 072	in (..) = Anzahl der Schneiden
	... F	Gelber Ring/Yellow ring/Bague jaune	Fein/Fine/Fin	ISO No. 041	in (..) = number of blades
	... UF	Weißer Ring/White ring/Bague blanche	Ultrafein/Ultrafine/Ultrafin	ISO No. 031	in (..) = nombre de lames